

Dimercaprol / BAL

Ph. Eur.

1. MANUFACTURER OF THE FINAL SUBSTANCE

Evonik Operations GmbH
Gutenbergstrasse 2
69221 Dossenheim
Germany

2. PRODUCT

- Nomenclature

Recommended International Non-proprietary name (INN):

Dimercaprol
British Anti-Lewisite (BAL)

Chemical name(s):


2,3-Dimercapto-1-propanol
2,3-Bis-sulfanyl-propan-1-ol (IUPAC)
1,2-Dithioglycerin
1-Propanol, 2,3-dimercapto- (CAS)
2,3-Dimercaptopropanol
2,3-Dimercaptopropan-1-ol
1,2-Thioglycerol

CAS No.: 59-52-9

Dimercaprol / BAL

Ph. Eur.

- Structure



Molecular formula: $C_3H_8OS_2$

Molecular mass: 124.2 g/mol

- General Properties

Description: Colorless or yellowish liquid with characteristic odor of mercaptans.

Solubility: Soluble in water (87 g/L) and peanut oil, miscible with benzyl benzoate.

Refractive index: $n_D^{20} = 1.568 - 1.574$

Dimercaprol / BAL

Ph. Eur.

3. SPECIFICATION

Table 1: Specification

Test Items	Test Methods	References	Acceptance Criteria	Frequencies
Appearance	Visual	Ph. Eur. monograph 0389	Clear, colorless to slightly yellow liquid	Each batch
Identity	De-colorization of iodine solution		conforms	
	Copper sulfate precipitation		conforms	
	Sodium Bismuthate-colorization		conforms	
Assay	Titration		98.5 - 101.5%	
Clarity	Turbidity	Ph. Eur. 2.2.1	< 3 NTU (clear)	
Color	Colorimetry	Ph. Eur. monograph 0389	not more intensely colored than reference solution B ₆ or BY ₆	
Acidity or alkalinity	---		nmt. 0.5 mL 0.01 N NaOH per 200 mg sample	
Refractive index	---		$n_D^{20} = 1.568 - 1.574$	
Halogenides (Halides)	Titration		≤ 1.0 mL 0.1 N NH ₄ SCN difference per 2 g sample	

Dimercaprol / BAL

Ph. Eur.

Table 1: Specification (cont.)

Test Items	Test Methods	References	Acceptance Criteria	Frequencies
Impurities:				
Impurity RRT 0.57	Gas Chromatography (GC-1)	Ph. Eur. 2.2.28	nmt. 0.5%	Each batch
Impurity RRT 0.60			nmt. 0.2%	
Impurity RRT 0.91			nmt. 0.5%	
Impurity 2,3-Dibromo-1-propanol			nmt. 0.5%	
Impurity 1,2,3-Trimercaptopropane			nmt. 1.5%	
Impurity RRT 1.43			nmt. 0.2%	
Impurity RRT 1.44			nmt. 0.2%	
Impurity RRT 1.82			nmt. 0.2%	
Impurity RRT 1.85			nmt. 0.2%	
Other impurities, each			nmt. 0.1%	
Others impurities, sum of all			nmt. 0.5%	
Residual solvents:				
Chloroform	Gas Chromatography (GC-2)	Ph. Eur. 2.2.28	nmt. 60 ppm	Each batch
Methanol			nmt. 3000 ppm	
Acetone			nmt. 5000 ppm	
Dimethyl formamide (DMF)	Gas Chromatography (GC-3)		nmt. 880 ppm	

- nmt. = not more than

Dimercaprol / BAL

Ph. Eur.

4. PACKAGING

- Container Closure System I

In narrow-neck bottles or jerry cans, depending on the amount and customer request. Generally, 5 L jerry cans or 0.1 L and 1 L narrow neck bottles are used.

Narrow-neck bottles: 0.1 L or 1 L

High Density Polyethylene (HDPE)

Jerry cans: 5 L

Inner layer: High Density Polyethylene (HDPE)

Outer layer: Ethylene – Vinyl alcohol co-polymer, EVOH (antistatic coating)

- Container Closure System II

Fiber drums

Compliance

The packaging is certified according to directives of EU and FDA, it conforms to the "Food Contact Regulations".

Material conforms to EU-framework on materials and articles intended for food contact: (EC) No. 1935/2004.

Raw material conforms to

- Regulation (EU) No. 10/2011 and amendments
- Commodities regulation of 10.04.1992 and amendments

All applied packaging systems have an "UN"-approval for transport of dangerous goods.

5. STORAGE AND STABILITY

Based on stability data, a re-test period of 36 months, if stored at 2 – 8 °C in the original HDPE containers, is defined.

Dimercaprol / BAL

Ph. Eur.

6. GENERAL INFORMATION

- GMP certificate: available.
- CEP No.: CEP 2007-330
- Residual solvents

The residual organic solvents are controlled in the final drug substance. The acceptance criteria are compliant with the guideline ICH Q3C.

Table 2: Residual solvents

Solvent Class 1	Acceptance criterion	Tested in routine	Method used
--	--	--	--
Solvent Class 2	Acceptance criterion	Tested in routine	Method used
Chloroform	nmt. 60 ppm	yes	GC
Methanol	nmt. 3000 ppm		
Dimethyl formamide	nmt. 880 ppm		
Solvent Class 3	Acceptance criterion	Tested in routine	Method used
Acetone	nmt. 5000 ppm	yes	GC

- nmt. = not more than

This information and any recommendations, technical or otherwise, are presented in good faith and believed to be correct as of the date prepared. Recipients of this information and recommendations must make their own determination as to its suitability for their purposes. In no event shall Evonik assume liability for damages or losses of any kind or nature that result from the use of or reliance upon this information and recommendations. Evonik expressly disclaims any representations and warranties of any kind, whether express or implied, as to the accuracy, completeness, non-infringement, merchantability and/or fitness for a particular purpose (even if Evonik is aware of such purpose) with respect to any information and recommendations provided. Reference to any trade names used by other companies is neither a recommendation nor an endorsement of the corresponding product, and does not imply that similar products could not be used. Evonik reserves the right to make any changes to the information and/or recommendations at any time, without prior or subsequent notice.