

# Determination of the pH value of hydrogen peroxide

## GENERAL INFORMATION ABOUT THE METHOD

This method describes the determination of the pH of hydrogen peroxide. The measurement is carried out directly in the present hydrogen peroxide sample. For the measurement, a commercial pH glass electrode and a commercial pH meter are required.

## EQUIPMENT

- pH meter
- pH glass electrode
- thermostat
- beakers, 50ml

## REAGENTS

- hydrogen peroxide solution (testing material)
- high purity water (osmosis and ion exchange treated drinking water)
- buffer solutions (e.g. pH 2, pH 4 and pH 7)
- storage solution for pH electrode (individually according to manufacturer's instructions)

## SPECIAL SAFETY INSTRUCTIONS

All reagents and chemicals must be handled according to the health and safety regulations. Refer to the safety data sheets.

## SPECIAL PROCEDURE INSTRUCTIONS

Danger of decomposition by contact with incompatible materials, contaminants, metals, alkalis, reducing agents.

## PROCEDURE

The equipment (pH meter and pH glass electrode) should be checked before each measurement. Therefore, perform a calibration with 2 or 3 different buffer solutions. The pH values of the buffer solutions must be known exactly. The selection of the buffer solutions depends on the expected pH of the hydrogen peroxide sample and of the pH meter type (observe manufacturer's calibration instructions!).

After the calibration temper the hydrogen peroxide sample to 20°C.

To measure the sample pH, remove the electrode from the storage tube and rinse it with high purity water. The glass membrane should not be touched, a possibly attached drop pick up carefully with a clean tissue.

Immerse the electrode in the sample solution and wait until the indicated pH value has stabilized. Then read off and note the value.

## CALCULATION

Not applicable

## ENVIRONMENT/DISPOSAL OF CHEMICALS

The disposal of laboratory quantities of hydrogen peroxide must be in accordance with local regulations.

## LITERATURE

- Manufacturers equipment descriptions
- Product information "Hydrogen Peroxide"

## REMARKS

The method is based on the internal analytical method WM30.

# Determination of the pH value of hydrogen peroxide

## Disclaimer

This information and all further technical advice are based on our present knowledge and experience. However, it implies no liability or other legal responsibility on our part, including with regard to existing third-party intellectual property rights, especially patent rights. In particular, no warranty, whether express or implied, or guarantee of product properties in the legal sense is intended or implied. We reserve the right to make any changes according to technological progress or further developments. The customer is not released from the obligation to conduct careful inspection and testing of incoming goods. Performance of the product described herein should be verified by testing, which should be carried out only by qualified experts in the sole responsibility of a customer. Reference to trade names used by other companies is neither a recommendation, nor does it imply that similar products could not be used.

Evonik Operations GmbH  
Active Oxygens  
Rodenbacher Chaussee 4  
63457 Hanau, Germany  
PHONE +49 6181 59-3295  
FAX +49 6181 59-73295  
[www.active-oxygens.com](http://www.active-oxygens.com)

**FUTURIZE PEROXIDE**  
The peroxide experts at Evonik